

# Analysis V3 shortcuts

## Mouse events

---

Context menu	Right click
Pan / move spectrum around	Left drag
Zoom into area	Shift right drag
Zoom into area	Middle drag
Zoom in/out around cursor	display
Zoom in/out for selected axis	Scroll wheel in axis active area
Select new	Left click
Add/Remove one to/from Selection	Ctrl/Cmd left click
Select in area and add	Ctrl/Cmd left drag
Pick one peak at position	Ctrl/Cmd shift left click
Pick peaks in area	Ctrl/Cmd shift left drag
Move single peak	Ctrl/Cmd middle drag
Move module	Drag&drop purple bar
Move whole window	Alt Left drag
Pop out module (manage as windows)	Double click on Purple bar
Make new object	in Sidebar
View/Edit object	in Sidebar

## Keyboard shortcuts

---

### *Graphics :*

	Shortcut	Related
Delete selected	Del (Fn+Backspace)	
Toggle Grid	GS	CH, TA, PC
Toggle Crosshairs	CH	GS, TA, PC
Toggle Double Crosshairs	CD	CH
Toggle Horizontal Trace	TH	TV, TU, TD
Toggle Vertical Trace	TV	TH, TU, TD
Increase Trace Scale	TU	TH, TV, TD
Decrease Trace Scale	TD	TH, TV, TU
Create Mark at cursor position	MK	MC
Clear All Marks	MC	MK
Navigate to Peak position	FP	FN
Navigate to NmrResidue position	FN	FP
Copy PeakList	CL	CP
Copy Peaks	CP	CL
Clear Current Selected Peaks	CZ	
Refit Selected Peaks	RP	SE
Snap to Extremum Current Selected Peaks	SE	RP
Select All Peaks	Ctrl/Cmd+A	
Set Phasing Pivot to cursor position	PV	PC, TA, TR
Add Integral in 1D plot	I1	
Add Phasing Trace	TA	PC, PV, TR
Remove Phasing Trace	TR	PC, PV, TA
Stack 1D spectra	SK	

Switch Mouse Mode	MM	
Zoom in to Spectrum	ZI or +	ZO, ZA
Zoom out of Spectrum	ZO or -	ZI, ZA
Zoom to Full Spectrum	ZA	ZI, ZO
Store Zoom setting	ZS	ZR, ZP, ZN, ZA
Restore stored Zoom setting	ZR	ZS, ZP, ZN, ZA
Go to previous Zoom level	ZP	ZN, ZS, ZR, ZA
Go to next Zoom level	ZN	ZP, ZS, ZR, ZA
Toggle Share Last Axis Across Display	LA	GS, CH, PC
Cycle Peak Labelling	PL	PS
Cycle Peak Symbols	PS	PL
Set Contour Levels	CO	

*Project :*

New project	Ctrl/Cmd+N	
Open project	Ctrl/Cmd+O	
Load Data	LD	
Save	Ctrl/Cmd+S	
Save as	SA	
Import Nef File	IN	
Export Nef File	EX	
Undo	Ctrl/Cmd+Z	
Redo	Ctrl/Cmd+Y	
Preferences	Ctrl/Cmd+,	
Quit	Ctrl/Cmd+Q	

*Spectrum :*

Load Spectrum	LS	SS
Show Spectrum Groups	SS	LS
Set Experiment Types	ET	
Pick 1D Peaks	P1	PP
Pick ND Peaks	PP	P1
Copy PeakList	CL	CP
Copy Peaks	CP	CL
Make Projection	PJ	
Print to File	Ctrl/Cmd+P	

*Molecules :*

Residue information	RI
Reference Chemical Shifts	RC

*Assign :*

Set up NmrResidues	SN
Pick and assign	PA
Backbone assignment	BB
Peak Assigner	AA
Assignment Inspector	AI

*View :*

NmrAtom Chemical Shift Table	CT
------------------------------	----

NmrResidue Table	NT	
Peak Table	PT	
Integral Table	IT	
Restraint Table	RT	
Structure Table	ST	
Sequence Graph	SG	
Atom Assigner	AS	
Chemical Shift Mapping	CM	
Notes Table	NO	
Python Console	Space-Space	
<i>For current strip</i>		
Show/Hide Toolbar	TB	SB
Show/Hide Spectrum Toolbar	SB	TB
Show/Hide Phasing Console	PC	PV, TA, TR
Move to Next Spectrum (and Hide Current)	TabTab	TabQ
Move to Previous Spectrum (and Hide Current)	TabQ	TabTab
Show all Spectra on Current Display	TabA	TabZ
Hide all Spectra on Current Display	TabZ	TabA
Copy with X-Y Axes Flipped	XY	XZ, YZ
Copy with X-Z Axes Flipped	XZ	XY, YZ
Copy with Y-Z Axes Flipped	YZ	XY, XZ
Estimate Noise	EN	
 <i>Macro :</i>		
Define User Shortcuts	DU	
Reserved for user macros	Q0, Q1, ...Q9	
Reserved for user macros	W0, W1, ...W9	
Reserved for user macros	Y0, Y1, ...Y9	