

Analysis V3 shortcuts

Mouse events

| | |
|------------------------------------|----------------------------------|
| Context menu | Right click |
| Pan | Left drag |
| Zoom into area | Shift right drag |
| Zoom into area | Middle drag |
| Zoom in/out around cursor | display |
| Zoom in/out for selected axis | Scroll wheel in axis active area |
| Select new | Left click |
| Add/Remove one to/from Selection | Ctrl/Cmd left click |
| Select in area and add | Ctrl/Cmd left drag |
| Pick one peak at position | Ctrl/Cmd shift left click |
| Pick peaks in area | Ctrl/Cmd shift left drag |
| Move single peak | Ctrl/Cmd middle drag |
| Move module | Drag&drop purple bar |
| Move whole window | Alt Left drag |
| Pop out module (manage as windows) | Double click on Purple bar |
| Make new object | in Sidebar |
| View/Edit object | in Sidebar |

Keyboard shortcuts

Graphics :

| | Shortcut | Related |
|---|--------------------|------------|
| Delete selected | Del (Fn+Backspace) | |
| Toggle Grid | GS | CH, TA, PC |
| Toggle Crosshairs | CH | GS, TA, PC |
| Toggle Double Crosshairs | CD | CH |
| Toggle Horizontal Trace | TH | TV, TU, TD |
| Toggle Vertical Trace | TV | TH, TU, TD |
| Increase Trace Scale | TU | TH, TV, TD |
| Decrease Trace Scale | TD | TH, TV, TU |
| Create Mark at cursor position | MK | MC |
| Clear All Marks | MC | MK |
| Navigate to Peak position | FP | FN |
| Navigate to NmrResidue position | FN | FP |
| Copy PeakList | CL | CP |
| Copy Peaks | CP | CL |
| Clear Current Selected Peaks | CZ | |
| Refit Selected Peaks | RP | SE |
| Snap to Extremum Current Selected Peaks | SE | RP |
| Select All Peaks | Ctrl/Cmd+A | |
| Set Phasing Pivot to cursor position | PV | PC, PT, TR |
| Add Integral in 1D plot | I1 | |
| Add Phasing Trace | PT | PC, PV, TR |
| Remove Phasing Trace | TR | PC, PV, PT |
| Switch Mouse Mode | MM | |

| | | |
|---------------------------------------|---------|----------------|
| Zoom in to Spectrum | ZI or + | ZO, ZS, ZR, ZA |
| Zoom out of Spectrum | ZO or - | ZI, ZS, ZR, ZA |
| Zoom to Full Spectrum | ZA | ZI, ZS, ZR, ZA |
| Toggle Share Last Axis Across Display | LA | GS, CH, PC |
| Cycle Peak Labelling | PL | PS |
| Cycle Peak Symbols | PS | PL |

Project :

| | |
|-----------------|------------|
| New project | Ctrl/Cmd+N |
| Open project | Ctrl/Cmd+O |
| Load Data | LD |
| Save | Ctrl/Cmd+S |
| Save as | SA |
| Import Nef File | IN |
| Export Nef File | EX |
| Undo | Ctrl/Cmd+Z |
| Redo | Ctrl/Cmd+Y |
| Preferences | Ctrl/Cmd+, |
| Quit | Ctrl/Cmd+Q |

Spectrum :

| | | |
|----------------------|------------|----|
| Load Spectrum | LS | SS |
| Show Spectrum Groups | SS | LS |
| Set Experiment Types | ET | |
| Pick 1D Peaks | P1 | PP |
| Pick ND Peaks | PP | P1 |
| Copy PeakList | CL | CP |
| Copy Peaks | CP | CL |
| Make Projection | PJ | |
| Print to File | Ctrl/Cmd+P | |

Molecules :

| | |
|---------------------------|----|
| Reference Chemical Shifts | RC |
|---------------------------|----|

Assign :

| | |
|----------------------|----|
| Set up NmrResidues | SN |
| Pick and assign | PA |
| Backbone assignment | BB |
| Peak Assigner | AA |
| Assignment Inspector | AI |
| Residue information | RI |

View :

| | |
|------------------------------|----|
| NmrAtom Chemical Shift Table | CT |
| NmrResidue Table | NT |
| Peak Table | PT |
| Integral Table | IT |
| Restraint Table | RT |
| Structure Table | ST |
| Sequence Graph | SG |

| | | |
|---|---------------|----------------|
| Atom Assigner | AS | |
| Chemical Shift Mapping | CM | |
| Notes Table | NO | |
| Python Console | Space-Space | |
| <i>For current strip</i> | | |
| Show/Hide Toolbar | TB | SB |
| Show/Hide Spectrum Toolbar | SB | TB |
| Show/Hide Phasing Console | PC | PV, PT, TR |
| Move to Next Spectrum (and Hide Current) | TabTab | |
| Move to Previous Spectrum (and Hide Current) | TabQ | |
| Toggle all Spectra on Current Display | TabA | |
| Store Zoom | ZS | ZI, ZO, ZR, ZA |
| Restore Zoom | ZR | ZI, ZO, ZS, ZA |
| Copy with X-Y Axes Flipped | XY | XZ, YZ |
| Copy with X-Z Axes Flipped | XZ | XY, YZ |
| Copy with Y-Z Axes Flipped | YZ | XY, XZ |
| <i>Macro :</i> | | |
| Define User Shortcuts | DU | |
| Reserved for user macros | Q0, Q1, ...Q9 | |
| Reserved for user macros | W0, W1, ...W9 | |
| Reserved for user macros | Y0, Y1, ...Y9 | |