

# Analysis V3 shortcuts

## Mouse events

---

Context menu	Right click
Pan / move spectrum around	Left drag
Zoom into area	Shift right drag
Zoom into area	Middle drag
Zoom in/out around cursor	Scroll wheel in display
Zoom in/out for selected axis	Scroll wheel in axis active area
Select new	Left click
Add/Remove one to/from Selection	Ctrl/Cmd left click
Select in area and add	Ctrl/Cmd left drag
Pick one peak at position	Ctrl/Cmd shift left click
Pick peaks in area	Ctrl/Cmd shift left drag
Move single peak	Ctrl/Cmd middle drag
Move module	Drag&drop purple bar
Move whole window	Alt Left drag
Pop out module (manage as windows)	Double click on Purple bar
Make new object	Double click <New> in Sidebar
View/Edit object	Double click object in Sidebar

## Keyboard shortcuts

---

	Shortcut	Related
<i>Graphics :</i>		
Delete selected	Del (Fn+Backspace)	
Toggle Grid	GS	CH, TA, PC
Toggle Crosshairs	CH	GS, TA, PC
Toggle Double Crosshairs	CD	CH
Toggle Horizontal Trace	TH	TV, TU, TD
Toggle Vertical Trace	TV	TH, TU, TD
Increase Trace Scale	TU	TH, TV, TD
Decrease Trace Scale	TD	TH, TV, TU
Create Mark at cursor position	MK	MX, MY, MZ, MW, MC
Create Mark at cursor X-axis position	MX	MK, MY, MZ, MW, MC
Create Mark at cursor Y-axis position	MY	MK, MX, MZ, MW, MC
Create Mark at cursor Z-axis position	MZ	MK, MX, MY, MW, MC
Create Mark at cursor W-axis position	MW	MK, MX, MY, MZ, MC
Clear All Marks	MC	MK, MX, MY, MZ, MW PM, PX, PY, PZ, PW UM, UX, UY, UZ, UW
Create Mark through selected Peaks	PM	PX, PY, PZ, PW, MC
Create X-axis Mark through selected Peaks	PX	PM, PY, PZ, PW, MC
Create Y-axis Mark through selected Peaks	PY	PM, PX, PZ, PW, MC
Create Z-axis Mark through selected Peaks	PZ	PM, PX, PY, PW, MC
Create W-axis Mark through selected Peaks	PW	PM, PX, PY, PZ, MC
Create Mark through selected Multiplets	UM	UX, UY, UZ, UW, MC
Create X-axis Mark through selected Multiplets	UX	UM, UY, UZ, UW, MC
Create Y-axis Mark through selected Multiplets	UY	UM, UX, UZ, UW, MC
Create Z-axis Mark through selected Multiplets	UZ	UM, UX, UY, UW, MC

Create W-axis Mark through selected Multiplets	UW	UM, UX, UY, UZ, MC
--	----	--------------------

Navigate to Peak position	FP	FN
Navigate to NmrResidue position	FN	FP
Copy PeakList	CL	CP
Copy Peaks	CP	CL
Clear Current Selected Peaks	CZ	
Refit Selected Peaks	RP	SE
Snap to Extremum Current Selected Peaks	SE	RP
Select All Peaks	Ctrl/Cmd+A	
Set Phasing Pivot to cursor position	PV	PC, TA, TR
Add Integral in 1D plot	I1	
Add Phasing Trace	TA	PC, PV, TR
Remove Phasing Trace	TR	PC, PV, TA
Stack 1D spectra	SK	
Switch Mouse Mode	MM	
Zoom in to Spectrum	ZI or +	ZO, ZA
Zoom out of Spectrum	ZO or -	ZI, ZA
Zoom to Full Spectrum	ZA	ZI, ZO
Store Zoom setting	ZS	ZR, ZP, ZN, ZA
Restore stored Zoom setting	ZR	ZS, ZP, ZN, ZA
Go to previous Zoom level	ZP	ZN, ZS, ZR, ZA
Go to next Zoom level	ZN	ZP, ZS, ZR, ZA
Toggle Share Last Axis Across Display	LA	GS, CH, PC
Cycle Peak Labelling	PL	PS
Cycle Peak Symbols	PS	PL
Set Contour Levels	CO	

*Project :*

New project	md+N
Open project	md+O
Load Data	LD
Save	md+S
Save as	SA
Import Nef File	IN
Export Nef File	EX
Undo	md+Z
Redo	md+Y
Preferences	md+,
Quit	md+Q

*Spectrum :*

Load Spectrum	LS	SS
Show Spectrum Groups	SS	LS
Set Experiment Types	ET	
Pick 1D Peaks	P1	PP
Pick ND Peaks	PP	P1
Copy PeakList	CL	CP
Copy Peaks	CP	CL
Estimate Volumes	EV	
Make Projection	PJ	
Print to File	md+P	

*Molecules :*

Residue information	RI
Reference Chemical Shifts	RC

*Assign :*

Set up NmrResidues	SN
Pick and assign	PA
Backbone assignment	BB
Peak Assigner	AP
NmrAtom Assigner	AN
Assignment Inspector	AI

*View :*

NmrAtom Chemical Shift Table	CT
NmrResidue Table	NT
Peak Table	PT
Integral Table	IT
Restraint Table	RT
Structure Table	ST
Sequence Graph	SG
Chemical Shift Mapping	CM
Notes Table	NO
Python Console	Space

*For current strip :*

Show/Hide Toolbar	TB	SB
Show/Hide Spectrum Toolbar	SB	TB
Show/Hide Phasing Console	PC	PV, TA, TR
Move to Next Spectrum (and Hide Current)	b	TabQ
Move to Previous Spectrum (and Hide Current)	TabQ	TabTab
Show all Spectra on Current Display	TabA	TabZ
Hide all Spectra on Current Display	TabZ	TabA
Copy with X-Y Axes Flipped	XY	XZ, YZ
Copy with X-Z Axes Flipped	XZ	XY, YZ
Copy with Y-Z Axes Flipped	YZ	XY, XZ
Estimate Noise	EN	
Move up a z-plane	ZZ	XX
Move down a z-plane	XX	ZZ
Raise contour level	QQ	WW
Lower contour level	WW	QQ

*Macro :*

Define User Shortcuts	DU
Reserved for user macros	Q1,
Reserved for user macros	W1,
Reserved for user macros	...Y9