

Analysis V3 shortcuts

Mouse events

Context menu	Right click
Pan / move spectrum around	Left drag
Zoom into area	Shift right drag
Zoom into area	Middle drag
Zoom in/out around cursor	Scroll wheel in display
Zoom in/out for selected axis	Scroll wheel in axis active area
Select new	Left click
Add/Remove one to/from Selection	Ctrl/Cmd left click
Select in area and add	Ctrl/Cmd left drag
Pick one peak at position	Ctrl/Cmd shift left click
Pick peaks in area	Ctrl/Cmd shift left drag
Move single peak	Ctrl/Cmd middle drag
Move module	Drag&drop purple bar
Move whole window	Alt Left drag
Pop out module (manage as windows)	Double click on Purple bar
Make new object	Double click <New> in Sidebar
View/Edit object	Double click object in Sidebar

Keyboard shortcuts

	Shortcut	Related
<i>Graphics :</i>		
Delete selected	Del (Fn+Backspace)	
Toggle Grid	GS	CH, TA, PC
Toggle Crosshairs	CH	GS, TA, PC
Toggle Double Crosshairs	CD	CH
Toggle Horizontal Trace	TH	TV, TU, TD
Toggle Vertical Trace	TV	TH, TU, TD
Increase Trace Scale	TU	TH, TV, TD
Decrease Trace Scale	TD	TH, TV, TU
Create Mark at cursor position	MK	MX, MY, MZ, MW, MC
Create Mark at cursor X-axis position	MX	MK, MY, MZ, MW, MC
Create Mark at cursor Y-axis position	MY	MK, MX, MZ, MW, MC
Create Mark at cursor Z-axis position	MZ	MK, MX, MY, MW, MC
Create Mark at cursor W-axis position	MW	MK, MX, MY, MZ, MC
Clear All Marks	MC	MK, MX, MY, MZ, MW PM, PX, PY, PZ, PW UM, UX, UY, UZ, UW
Create Mark through selected Peaks	PM	PX, PY, PZ, PW, MC
Create X-axis Mark through selected Peaks	PX	PM, PY, PZ, PW, MC
Create Y-axis Mark through selected Peaks	PY	PM, PX, PZ, PW, MC
Create Z-axis Mark through selected Peaks	PZ	PM, PX, PY, PW, MC
Create W-axis Mark through selected Peaks	PW	PM, PX, PY, PZ, MC
Create Mark through selected Multiplets	UM	UX, UY, UZ, UW, MC
Create X-axis Mark through selected Multiplets	UX	UM, UY, UZ, UW, MC
Create Y-axis Mark through selected Multiplets	UY	UM, UX, UZ, UW, MC
Create Z-axis Mark through selected Multiplets	UZ	UM, UX, UY, UW, MC

Create W-axis Mark through selected Multiplets	UW	UM, UX, UY, UZ, MC
--	----	--------------------

Navigate to Peak position	FP	FN
Navigate to NmrResidue position	FN	FP
Copy PeakList	CL	CP
Copy Peaks	CP	CL
Clear Current Selected Peaks	CZ	
Refit Selected Peaks	RP	SE
Snap to Extremum Current Selected Peaks	SE	RP
Select All Peaks	Ctrl/Cmd+A	
Set Phasing Pivot to cursor position	PV	PC, TA, TR
Add Phasing Trace	TA	PC, PV, TR
Remove Phasing Trace	TR	PC, PV, TA
Add Integral in 1D plot	I1	
Stack 1D spectra	SK	
Switch Mouse Mode	MM	
Zoom in to Spectrum	ZI or +	ZO, ZA
Zoom out of Spectrum	ZO or -	ZI, ZA
Zoom to Full Spectrum	ZA	ZI, ZO
Store Zoom setting	ZS	ZR, ZP, ZN, ZA
Restore stored Zoom setting	ZR	ZS, ZP, ZN, ZA
Go to previous Zoom level	ZP	ZN, ZS, ZR, ZA
Go to next Zoom level	ZN	ZP, ZS, ZR, ZA
Set Zoom...	SZ	
Toggle Share Last Axis Across Display	LA	GS, CH, PC
Cycle Peak Labelling	PL	PS
Cycle Peak Symbols	PS	PL
Set Contour Levels	CO	

File :

New project	Ctrl/Cmd+N
Open project	Ctrl/Cmd+O
Load Data	LD
Save	Ctrl/Cmd+S
Save as	SA
Import Nef File	IN
Import BMRB File	BI
Export Nef File	EX
Preferences	Ctrl/Cmd+,
Quit	Ctrl/Cmd+Q

Edit :

Undo	Ctrl/Cmd+Z
Redo	Ctrl/Cmd+Y

View :

NmrAtom Chemical Shift Table	CT
NmrResidue Table	NT
Peak Table	PT
Integral Table	IT
Multiplet Table	MT
Restraint Table	RT
Structure Table	ST

Data Table	DT
Sequence Graph	SG
Violation Table	VT
Restraint Analysis Table	AT
Chemical Shift Mapping	CM
Notes Editor	NO
Python Console	Space-Space

In Active Spectrum Display / Strip :

Show/Hide Toolbar	TB	SB
Show/Hide Spectrum Toolbar	SB	TB
Show/Hide Phasing Console	PC	PV, TA, TR
Set Zoom	IT	
Restore Zoom	MT	
Move to Next Spectrum (and Hide Current)	TabTab	TabQ
Move to Previous Spectrum (and Hide Current)	TabQ	TabTab
Show all Spectra on Current Display	TabA	TabZ
Hide all Spectra on Current Display	TabZ	TabA
Copy with X-Y Axes Flipped	XY	XZ, YZ
Copy with X-Z Axes Flipped	XZ	XY, YZ
Copy with Y-Z Axes Flipped	YZ	XY, XZ
Estimate Noise	EN	
Move up a z-plane	JJ	KK
Move down a z-plane	KK	JJ
Raise contour level	QQ	WW
Lower contour level	WW	QQ

Spectrum :

Load Spectrum	LS	
Show Spectrum Groups	SS	
Validate Paths	VP	
Set Experiment Types	ET	
Pick 1D Peaks	P1	PP
Pick ND Peaks	PP	P1
Copy PeakList	CL	CP
Copy Peaks	CP	CL
Estimate Volumes	EV	
Estimate Current Peak Volumes	EC	
Reorder Peak List Axes	RL	
Make Strip Plot	SP	
Make Projection	PJ	
Print to File	Ctrl/Cmd+P	

Molecules :

Residue information	RI
Reference Chemical Shifts	RC

Assign :

Set up NmrResidues	SN
Pick and assign	PA
Backbone assignment	BB
Peak Assigner	AP

NmrAtom Assigner
Assignment Inspector

AN
AI

Macro :

New Macro

NM OM, RM

Open User Macro

OM NM, RM

Run Macro

RM NM, OM

Define User Shortcuts

DU

Reserved for user macros

Q0, Q1, ...Q9

Reserved for user macros

W0, W1, ...W9

Reserved for user macros

Y0, Y1, ...Y9