

Analysis V3_{beta2} shortcuts

Mouse events

Context menu	Right click
Pan	Left drag
Zoom into area	Shift right drag
Zoom into area	Middle drag
Zoom in/out around cursor	Scroll wheel in display
Zoom in/out for selected axis	Scroll wheel in axis active area
Select new	Left click
Add/Remove one to/from Selection	Ctrl/Cmd left click
Select in area and add	Ctrl/Cmd left drag
Pick one peak at position	Ctrl/Cmd shift left click
Pick peaks in area	Ctrl/Cmd shift left drag
Move single peak	Ctrl/Cmd middle drag
Move module	Drag&drop purple bar
Move whole window	Alt Left drag
Pop out module (manage as windows)	Double click on Purple bar
Make new object	Double click <New> in Sidebar
View/Edit object	Double click object in Sidebar

Keyboard shortcuts

Graphics :

	Shortcut
Delete selected	Del (Fn+Backspace)
Toggle Grid	GS
Toggle Crosshairs	CH
Toggle Horizontal Trace	TH
Toggle Vertical Trace	TV
Increase Trace Scale	TU
Decrease Trace Scale	TD
Create Mark at cursor position	MK
Clear All Marks	MC
Navigate to Peak position	FP
Navigate to NmrResidue position	FN
Copy PeakList	CL
Copy Peaks	CP
Clear Current Selected Peaks	CZ
Refit Selected Peaks	RP
Snap to Extremum Current Selected Peaks	SE
Select All Peaks	Ctrl/Cmd+A
Toggle Phasing Pivot	PV
Add Integral in 1D plot	I1
Add Phasing Trace	PT

Remove Phasing Trace	TR
Switch Mouse Mode	MM
Zoom in to Spectrum	ZI
Zoom out of Spectrum	ZO
Toggle Share Y Axis Across Display	TA

Project :

New project	Ctrl/Cmd+N
Open project	Ctrl/Cmd+O
Load Data	LD
Save	Ctrl/Cmd+S
Save as	SA
Import Nef File	IN
Export Nef File	EX
Undo	Ctrl/Cmd+Z
Redo	Ctrl/Cmd+Y
Quit	Ctrl/Cmd+Q

Spectrum :

Load Spectrum	LS
Show Spectrum Groups	SS
Set Experiment Types	ET
Pick 1D Peaks	P1
Pick ND Peaks	PP
Copy PeakList	CL
Copy Peaks	CP
<i>Make Projection</i>	PJ
Print to File	PR

Molecules :

Reference Chemical Shifts	RC
---------------------------	----

Assign :

Set up NmrResidues	SN
Pick and assign	PA
Backbone assignment	BB
Peak Assigner	AA
Assignment Inspector	AI
Residue information	RI

View :

NmrAtom Chemical Shift Table	CT
NmrResidue Table	NT
Peak Table	LT
Integral Table	IT
Restraint Table	RT
Structure Table	ST

Sequence Graph	SG
Atom Selector	AS
Chemical Shift Mapping	CM
Show Sequence	SQ
Notes Table	NO
Python Console	PY
<i>For current strip</i>	
Show/Hide Toolbar	TB
Show/Hide Phasing Console	PC
Move to Next Spectrum (and Hide Current)	Double Tab
Move to Previous Spectrum (and Hide Current)	Double Backspace
Store Zoom	ZS
Restore Zoom	ZR
Flip X-Y Axis	XY
Flip X-Z Axis	XZ
Flip Y-Z Axis	YZ
 <i>Macro :</i>	
Define User Shortcuts	DU
Reserved for user macros	Q0, Q1, ...Q9
Reserved for user macros	W0, W1, ...W9
Reserved for user macros	Y0, Y1, ...Y9