

Analysis V3 shortcuts

Mouse events

	Shortcut
Context menu	Right click
Pan	Left drag
Zoom into area	Shift right drag
Zoom into area	Middle drag
Zoom in/out around cursor	Scroll wheel in display
Zoom in/out for selected axis	Scroll wheel in axis active area
Select new	Left click
Add/Remove one to/from Selection	Ctrl/Cmd left click
Select in area and add	Ctrl/Cmd left drag
Pick one peak at position	Ctrl/Cmd shift left click
Pick peaks in area	Ctrl/Cmd shift left drag
Move single peak	Ctrl/Cmd middle drag
Move module	Drag&drop purple bar
Move whole window	Alt Left drag
Pop out module (manage as windows)	Double click on Purple bar
Make new object	Double click <New> in Sidebar
View/Edit object	Double click object in Sidebar

Keyboard shortcuts

	Shortcut
<i>Graphics :</i>	
Delete selected	Del (Fn+Backspace)
Toggle Grid	GS
Toggle Crosshairs	CH
Toggle Horizontal Trace	TH
Toggle Vertical Trace	TV
Increase Trace Scale	TU
Decrease Trace Scale	TD
Create Mark at cursor position	MK
Clear All Marks	MC
Navigate to Peak position	FP
Navigate to NmrResidue position	FN
Copy PeakList	CL
Copy Peaks	CP
Clear Current Selected Peaks	CZ
Refit Selected Peaks	RP
Snap to Extremum Current Selected Peaks	SE
Select All Peaks	Ctrl/Cmd+A

Project :

New project	PN
Open project	PO
Load Spectrum	LS
Load Data	LD
Save	PS
Save as	SA
Undo	Ctrl/Cmd+Z
Redo	Ctrl/Cmd+Y
Quit	QT

Spectrum :

<i>Spectrum Groups</i>	SS
Set Experiment Types	ET
Pick Peaks	PP
<i>Peak Integration (1D)</i>	IT
<i>Make Projection</i>	PJ
Print to File	PR

Molecules :

Create new molecule	CM
Show sequence	SQ
Reference Chemical Shifts	RC

Assign :

Set up NmrResidues	SN
Pick and assign	PA
Backbone assignment	BB
Peak Assigner	AA
Assignment Inspector	AI
Residue information	RI

View :

New Blank Display	ND
NmrAtom Chemical Shift Table	CT
NmrResidue Table	NT
Peak Table	LT
Sequence Graph	SG
Atom Selector	AS
Python Console	PY
<i>For current strip</i>	
Show/Hide Toolbar	TB
Show/Hide Phasing Console	PC
Move to Next Spectrum (and Hide Current)	MN
Move to Previous Spectrum (and Hide Current)	MP
Reset Zoom	RZ

Macro :

Run macro

Reserved for user macros

Reserved for user macros

Reserved for user macros

RM

Q1, Q2, ...Q0

W1, W2, ...W0

U1, U2, ...U0