

Analysis V3 shortcuts

Mouse events

| Action | Shortcut |
|------------------------------------|----------------------------------|
| Raise Context menu | Right click |
| Move modules | Drag&drop purple bar |
| Move whole window | Alt Left click and drag |
| Pan | Left click and drag |
| Pick peaks | Ctrl left drag |
| Pop out module (manage as windows) | Double click on Purple bar |
| Select | Left click |
| Multi-Select | Ctrl Left click |
| Select area | Shift left click and drag |
| Spectrum information | Double click ('SP' in Sidebar) |
| Show PeakList | Double click ('PL' in Sidebar) |
| Zoom in/out around cursor | Scroll wheel in display |
| Zoom in/out for selected axis | Scroll wheel in axis active area |

Keyboard shortcuts

| Action | Shortcut |
|--------------------------------|------------------|
| New project | PN |
| Open project | PO |
| Load Data | LD |
| Save | PS |
| Save as | SA |
| Undo | Ctrl+Z |
| Redo | Ctrl+Y |
| Close program | QT |
| Peak Table | LT |
| Find Peaks | PP |
| Delete peaks | Delete (fn+Bksp) |
| Load Sequence | LS |
| Show sequence | SQ |
| Pick and assign | PA |
| Backbone assignment | BB |
| Residue information | RI |
| New Blank Display | ND |
| Run macro | RM |
| Show Console | PY |
| Show Reference Chemical Shifts | RC |

| | |
|---------------------------------|----|
| Show Assignment Module | AA |
| Show Atom Selector | AS |
| Show NMR Residue Table | NR |
| Toggle Grid | GS |
| Toggle Crosshairs | CH |
| Create Mark at cursor position | MK |
| Clear All Marks | MC |
| Navigate to Peak position | FP |
| Navigate to NmrResidue position | FN |
| Copy Assignments | CA |
| Show Experiment Type Popup | ET |
| Increase Trace Scale | TU |
| Decrease Trace Scale | TD |
| Show Horizontal Trace | TH |
| Show Vertical Trace | TV |
| Toggle Phasing Console | PC |
| Add Phasing Traces | PT |
| Remove Phasing Traces | PR |